### Modulus of elasticity



Physics	Mechanics	Fabric &	material properties
Applied Science	Engineering	Applied Mechanics	
Applied Science	Engineering	Materials Science	Mechanical Properties
Applied Science	Medicine	Biomechanics	
C Difficulty level	<b>R</b> Group size	C Preparation time	Execution time
easy	2	45+ minutes	45+ minutes







# **General information**

### **Application**







Fig.1: Experimental set-up

Solid matter changes form under the application of a force. The exact knowlege about the behaviour of such changes is very important in construction. The modulus of elasticity is used to describe this behaviour in mathematical way.







www.phywe.de

### **Theory (1/4)**

**DHIME** 

excellence in science

If a body is regarded as a continuum, and  $\vec{r}_0$  or  $\vec{r}$  defines the local vector of a point P in its undeformed and deformed state respectively, then, for small displacement vectors

$$ec{u}=ec{r}-ec{r}_0\equiv (u_{1,}u_{2,}u_{3})$$
 .

the deformation tensor  $\vec{d}$  is:

$$d_{\mathrm{ik}} = rac{\partial u_i}{\partial x_k} + rac{\partial u_k}{\partial x_i}$$

The forces  $d\vec{F}$  which act upon an element of volume of the body whose edges are intersected parallel to the coordinate surfaces will be described by the stress tensor  $\hat{\tau}$ 

This allocates a stress  $\vec{p}$  to each area element dA defined by the unit vector  $\vec{e}$  in the direction of the normal.

$$\overrightarrow{p} = rac{d\overrightarrow{F}}{dA}$$
 $\overrightarrow{p} = \overrightarrow{e} \odot \overrightarrow{ au}$ 

From Hooke's Law we optain the relationship between  $\hat{d}$  and  $\hat{\tau}$ :

$$au_{
m ik} = \sum_{
m I,m} c_{
m ik}^{l,m} d_{
m lm}$$

### **Theory (2/4)**

The tensor  $\hat{e}$  is symmetrical for an elastic body so that of the 81 components only 21 remain. This number is reduced to 2 for the isotropic elastic body, i.e. the modulus of elasticity E and either the shear modulus G or the Poisson ratio  $\mu$ :

$$au_{11} = rac{E}{1+\mu} (d_{11} + rac{\mu}{1-2\mu} (d_{11} + d_{22} + d_{33}))$$

$$au_{12} = Gd_{12} = rac{1}{2} \cdot rac{E}{1+\mu} \cdot d_{12}$$

and analogously for  $\tau_{22}, \tau_{33}, \tau_{13}, \tau_{23}$ 

If a force acts in only one direction then

$$au_{22}= au_{33}=0$$

so we obtain (2)

$$\tau_{11} = E \cdot d_{11}$$

If a bar of height b and width a, supported at both ends by supports (separated by a distance L), is subjected to a force  $F_y$  acting at its centre, it behaves like a bar supported in the middle, its two ends being subjected to a force  $F_y/2$  in the opposing direction. In order to express the bending  $\lambda$  as a function of the modulus of elasticity E, let us first consider an element of volume



### **Theory (3/4)**



PHYWE excellence in science


 $\mathrm{d}V = \mathrm{d} \cdot a \cdot b$ 


the upper layer of which is shortened on bending, and the bottom layer lengthened. The length of the central layer remains unchanged (neutral fibre).

In fig. 2, I and II denote the sides before and after deformation.

Using the symbols given in fig. 2, we obtain:

 $d\lambda = x \cdot d\varphi = \frac{2\sigma x}{h}$ 





### Theory (4/4)

The elastic force  $dF_x$  which produces the extension dl, according to (1), is


$$\frac{\mathrm{dF}_x}{\mathrm{ds}} = E \frac{\mathrm{dl}}{\mathrm{dx}}$$

where  $ds = a \cdot dy$  is the area of the rotated layer.

The force produces a torque

$$\mathrm{d}T_z = y\mathrm{d}F_x = rac{2Ea\sigma}{b\cdot\mathrm{d}x}y^2\mathrm{d}y$$

The sum of these torques produced by the elastic forces must be equal to the torque produced by the external force  $F_y/2$ .



from which we obtain

$$\mathrm{d}\lambda = rac{6F_yx^2}{\mathrm{Eab}^3}\mathrm{d}x$$

and, after integration, the total bending

$$\lambda = rac{1}{4} \cdot (rac{L}{B})^3 \cdot rac{1}{a} \cdot rac{F_y}{E}.$$



#### Equipment

Position	Material	Item No.	Quantity
1	Dial gauge 10/0.01 mm	03013-00	1
2	Holder for dial gauge	03013-01	1
3	Flat bars, set	17570-00	1
4	Knife-edge with stirrup	03015-00	1
5	Bolt with knife-edge	02049-00	2
6	Weight holder, 10 g	02204-00	1
7	Spring balance,transparent, 1 N	03065-02	1
8	Tripod base PHYWE	02002-55	2
9	Support rod, stainless steel, I = 250 mm, d = 10 mm	02031-00	2
10	Support rod, stainless steel, 750 mm	02033-00	1
11	Right angle clamp expert	02054-00	5
12	Slotted weight, black, 10 g	02205-01	10
13	Slotted weight, black, 50 g	02206-01	6
14	Vernier calliper stainless steel 0-160 mm, 1/20	03010-00	1
15	Measuring tape, I = 2 m	09936-00	1
16	Fish line, I. 100m	02090-00	1





## **Setup and Procedure**

#### Setup

The set up is as shown in Fig. 1. The dial gauge is mounted on the knife edge with stirrup. The flat bars must be positioned accurately on the two supporting knife edges having the possibility to move in x and y direction. The geometric data of the set up and the bars must be recorded several times or at different positions.

As the dial gauge possesses a restoring force which obeys Hooke's Law F  $\sim$  s, its characteristic curve must first of all be determined.





### Procedure (1/2)

PHYWE excellence in science

PHYWE excellence in science

The resulting force Fr of the dial gauge is the sum of

 $F_r = F_h + F_f$  where

 $F_h$  = Static friction force (constant)

 $F_f$  = restoring force (F  $\sim$  s)

Since the static friction force always acts opposed to the direction of motion, care should be taken of constant direction of force while recording of the characteristic curve of the dial gauge as well as during proformance of the experiment.

Additionally, the plunger is manually raised (the probe is unloaded) and then gently lowered. Due to this procedure the resulting force F is:

 $F_r = F_f - F_h$ 

#### Characteristic curve of the dial gauge:

Dynamometer and plunger of the dial gauge are so mounted, that the dial gauge shows full deflection.

Through a reduction of the tension of force tof the dynamometer a force results according to the above conditions.

### Procedure (2/2)

With the help of the spring balance the characteristic curve of the dial gauge has to be recorded.

During experimentation the forces must be corrected accordingly.

This means the resulting force F must be calculated.

Hence, the effective force is the sum of the weights of the additional masses and the resulting force of the dial gauge.




Fig. 3: Characteristic curve of a dial gauge.



### **Evaluation**

### Results (1/3)

Table 1: The modulus of elasticity for different materials
--

Material	Dimensions [mm]	$E[N \cdot m^{-2}]$
Steel	10×1.5	2.059·10 <sup>11</sup>
Steel	10×2	2.063·10 <sup>11</sup>
Steel	10×3	2.171·10 <sup>11</sup>
Steel	15×1.5	2.204·10 <sup>11</sup>
Steel	20×1.5	2.111·10 <sup>11</sup>
Aluminium	10×2	6.702·10 <sup>10</sup>
Brass	10×2	9.222·10 <sup>10</sup>





Fig. 4: Bending of a bar as a function of the force (steel, L = 0.48 m, a = 10 mm, b = 1.5 mm).



PHYWE excellence in science



### Results (3/3)




Fig. 7: Bending of a bar as a function of its thickness, at constant force (steel, width = 10 mm).



Robert-Bosch-Breite 10 37079 Göttingen Tel.: 0551 604 - 0 Fax: 0551 604 - 107 PHYWE excellence in science