

Surface tension of liquids



Phy

Physics

Mechanics

Mechanics of liquids & gases



Difficulty level

medium



Group size

2



Preparation time

45+ minutes




Execution time

45+ minutes



General information

Application



Water dripping from a tap

Several effects of surface tension can be seen in our daily life, for examples:

- **Water droplets:** droplets of water tend to be pulled into a spherical shape by the cohesive forces of the surface layer.
- **Surface tension disinfectants:** Low surface tension allows them to spread out on the cell walls of bacteria and disrupt them.
- **Soaps and detergents:** Cleaning of clothes by lowering the surface tension of the water so that it more readily soaks into pores and soiled areas.

Other information (1/2)

PHYWE
excellence in science

Prior knowledge



Cohesion is the property of like molecules sticking together, being mutually attractive. Adhesion is the tendency of dissimilar particles or surfaces to cling to one another. For example, at liquid-air interfaces, surface tension results from the greater attraction of liquid molecules to each other (due to cohesion) than to the molecules in the air (due to adhesion).

Scientific principle



The cohesive forces in a liquid generate tension on its surface, the so-called surface tension. A metal ring that is plunged into a liquid is withdrawn from the liquid. At a certain tensile force, the liquid film will be disrupted from the ring. Based on the tensile force and ring diameter, the surface tension of a liquid can be calculated.

Other information (2/2)

PHYWE
excellence in science

Learning objective



Understanding the surface tension of a liquid.

Tasks



Examine the surface tension of a liquid.

Safety instructions

PHYWE
excellence in science

For this experiment the general instructions for safe experimentation in science lessons apply.

For H- and P-phrases please consult the safety data sheet of the respective chemical.

Theory (1/2)

PHYWE
excellence in science

A metal ring that is plunged into a liquid is withdrawn from the liquid. At a certain tensile force, the liquid film will be disrupted from the ring. Based on the tensile force and ring diameter, the surface tension of a liquid can be calculated by way of the following formula:

$$\rho = \frac{F - F_G}{2\pi d} = \text{surface tension in N/m} \quad (1)$$

F = indicated force in N immediately prior to the disruption

F_G = weight of the ring in N

d = diameter of the ring in m

Theory (2/2)

Note

The surface tension is always specific to the liquid that is used. In order to demonstrate this, the experiment can be performed with several different liquids. However, liquids with a very low surface tension should not be used, since otherwise the rather simple measurement set-up would lead to excessively large measurement errors.

Equipment

| Position | Material | Item No. | Quantity |
|----------|--------------------------------------|----------|----------|
| 1 | Tripod base PHYWE | 02002-55 | 1 |
| 2 | Right angle clamp expert | 02054-00 | 1 |
| 3 | Support rod, stainless steel, 750 mm | 02033-00 | 1 |
| 4 | Rod with hook | 02051-00 | 1 |
| 5 | Surface tension measuring ring | 17547-00 | 1 |
| 6 | Petri dish, d 150 mm, glass | 64757-00 | 1 |
| 7 | Water, distilled 5 l | 31246-81 | 1 |
| 8 | Lab jack, 200 x 200 mm | 02074-01 | 1 |
| 9 | Spring balance, transparent, 0,2N | 03065-01 | 1 |

PHYWE
excellence in science

Setup and procedure

Setup (1/2)

PHYWE
excellence in science

Experimental setup

The picture shows the experiment set-up.

Procedure

PHYWE
excellence in science

The force gauge is suspended on a rod with a hook. By way of a right-angle clamp, the rod is fastened to a support rod that is secured in the tripod base. The force gauge is adjusted to zero. Then, the measurement ring is suspended in the eyelet of the force gauge.

The cleaned Petri dish is filled with pure water and placed on the adjustable table. The height of the table is adjusted so that the ring plunges approximately 3 mm deep into the liquid.

Then, the adjustable table is slowly lowered until the liquid film disrupts from the ring. The maximum value of the force F immediately prior to the disruption of the liquid film is read off the force sensor. Then, the weight of the ring in air is determined.

Evaluation (1/3)

PHYWE
excellence in science

The surfaces of liquids have a certain surface tension due to the cohesive forces of the liquids. The diameter of the ring is determined and the surface tension is calculated based on equation (1). The following values represent a measurement example:

Tensile Force: $F = 0.080 \text{ N}$

Weight of the ring: $F_G = 0.056 \text{ N}$

Diameter of the ring: $d = 0.052 \text{ m}$

$$\rho = \frac{F - F_G}{2\pi d} = \frac{0.080 - 0.056}{2 \cdot 0.052 \cdot \pi} \frac{\text{N}}{\text{m}} = 0.073 \frac{\text{N}}{\text{m}}$$

The surface tension of water is approximately 0.073 N/m .

Evaluation (2/3)

Which of the following contribute to surface tension?

- only cohesive forces
- both cohesive forces and adhesive forces
- only adhesive forces

✓ Check

Evaluation (3/3)

Fill in the blank:

Due to forces, molecules on the surface of a liquid are pulled downward to surrounding liquid molecules, while forces pull them upwards due to surrounding gaseous molecules acting. This force on the surface molecules is called .

adhesive

surface tension


cohesive

unbalanced

✓ Check

| Slide | Score/Total |
|---|-------------|
| Slide 13: Contribution to surface tension | 0/1 |
| Slide 14: Surface tension forces | 0/4 |

Total Score  0/5

 Show solutions

 Retry