


Magnetic field outside a straight conductor



The goal of this experiment is to investigate the magnetic field around a straight conductor.

Physics

Electricity & Magnetism

Magnetism & magnetic field



Difficulty level

hard



Group size

2



Preparation time

10 minutes




Execution time

10 minutes

PHYWE
excellence in science

General information

Application

PHYWE
excellence in science

Setup

Moving electric charges result in a magnetic field. The magnetic field around a straight conductor is the fundamental form, that such a magnetic field takes.

This experiment tries to determine the form of such a magnetic field.

Other information (1/2)

PHYWE
excellence in science

Prior knowledge



No prior knowledge is required.

Scientific principle



A current which flows through one or two neighbouring straight conductors produces a magnetic field around them. The dependences of these magnetic fields on the distance from the conductor and on the current are determined.

Other information (2/2)

PHYWE
excellence in science

Learning objective



The goal of this experiment is to investigate the magnetic field around a straight conductor.

Tasks



Determination of the magnetic field

1. of a straight conductor as a function of the current,
2. of a straight conductor as a function of the distance from the conductor,
3. of two parallel conductors, in which the current is flowing in the same direction, as a function of the distance from one conductor on the line joining the two conductors,
4. of two parallel conductors, in which the current is flowing in opposite directions, as a function of the distance from one conductor on the line joining the two conductors.

Theory (1/3)

Maxwell's 1st equation for the case when electric fields \vec{E} , variable with time, are absent,

$$\int_C \vec{B} \cdot d\vec{s} = \mu_0 \int_A \vec{j} \cdot d\vec{A} \quad (1)$$

together with Maxwell's 4th equation,

$$\int_A' \vec{B} \cdot d\vec{A} = 0 \quad (2)$$

provides the relationship between steady electric current I flowing through the area A ,

$$I = \int_A \vec{j} \cdot d\vec{A}$$

and the magnetic field \vec{B} which it produces.

Theory (2/3)

C is the boundary of A , A' is the enclosed area, \vec{j} is the electrical current density, μ_0

is the magnetic field constant, $\mu_0 = 1.26 \cdot 10^{-6} \text{ Vs/Am}$.

From (1) and (2) one obtains for a long straight conductor

$$|\vec{B}| = \frac{\mu_0}{2\pi} \cdot \frac{I}{|\vec{r}|} \quad (3)$$

where \vec{r} is the distance of the conductor from the point at which the magnetic field is measured.




Fig. 1: Contribution of a conductor section dl to the magnetic field at point Q .

Theory (3/3)

The direction of \vec{B} is \perp both to \vec{r} and to \vec{j} .

For a finite conductor one obtains, with the notation of Fig. 1:

$$d\vec{B} = \frac{1}{4\pi} \mu_0 \frac{I}{r^3} d\vec{l} \times \vec{r} \quad (\text{Biot-Savart})$$

and from this

$$|\vec{B}| = \frac{\mu I}{4\pi r} (\cos \phi_1 - \cos \phi_2)$$

Equipment

Position	Material	Item No.	Quantity
1	Current conductors, set of 4	06400-00	1
2	Coil, 6 turns	06510-00	1
3	Coil, 140 turns, 6 tapings	06526-01	1
4	Clamping device for iron cores	06506-00	1
5	Iron core, I-shaped, laminated	06500-00	1
6	Iron core, U-shaped, laminated	06501-00	1
7	PHYWE power supply, variable DC: 12 V, 5 A / AC: 15 V, 5 A	13540-93	1
8	PHYWE Teslameter, digital	13610-93	1
9	Hall probe, axial	13610-01	1
10	Current transformer/Clamp Ammeter adaptor	07091-10	1
11	Digital multimeter, 600V AC/DC, 10A AC/DC, 20 M Ω , 200 μ F, 20 kHz, -20°C... 760°C	07122-00	1
12	Meter scale, l = 1000 mm	03001-00	1
13	Barrel base expert	02004-00	1
14	Support rod, stainless steel, 500 mm	02032-00	1
15	Right angle clamp expert	02054-00	1
16	G-clamp	02014-00	2
17	Connecting cord, 32 A, 500 mm, yellow	07361-02	2
18	Universal clamp	37715-01	1



Setup and Procedure

Setup and Procedure (1/2)

The experimental set-up is arranged as shown in Fig. 2. The current transformer is used to measure the secondary current (20A... 120A). Since the primary and secondary current have a linear relationship, the primary current can also be measured. However, a calibration curve for primary/secondary current should then be recorded for each conductor. Because of the heating of the conductors, the current must be readjusted or a "warm-up time" must be allowed to elapse.




Fig. 2: Experimental setup

Setup and Procedure (2/2)

PHYWE
excellence in science

A phase displacement can occur between the “construction-kit” transformer and the magnetic field meter, giving the illusion of a “negative” magnetic field (minimum of the magnetic field indicator with increasing current). This can be eliminated by reversing the polarity of the primary of the transformer.

Higher short-time secondary currents can be achieved by connecting the constant and variable voltage in series on the power unit. Attention should be paid to the correct phase angle.

PHYWE
excellence in science



Evaluation

Results (1/5)

PHYWE
excellence in science

From the regression line to the measured values of Fig. 3 with the exponential statement

$$Y = A \cdot X^B$$

the exponent $B = 0.97 \pm 0.01$

and the slope $A = 52.91 \pm 0.01 \text{ A/mT}$

with (3) this gives $\mu_0 = 1.3 \cdot 10^{-6}$

Because of the small zero-deflection due to the instrument and the effect of the other conductor and the "construction kit" transformer, it is appropriate to carry out the measurement with small distances (up to approx. 3 cm) and with large currents (approx. 100 A).




Fig. 3: Relationship between current value and magnetic field of a long conductor (distance between conductor and measuring point: 1.1 cm).

Results (2/5)

PHYWE
excellence in science




Fig. 4: Magnetic field of a long conductor as a function of distance ($I = 100 \text{ A}$).




Fig. 5: Magnetic field of two conductors 1 and 2.

Results (3/5)

For the case of two parallel conductors in the z-direction, both carrying the same current I in the same direction ($p = 1$) or in opposite directions ($p = -1$), the superposition of the magnetic fields gives the components B_x and B_y of the magnetic field at point Q with the notation of Fig. 5.

$$B_x = |\vec{B}_1| \sin \alpha_1 + p \cdot |\vec{B}_2| \sin \alpha_2 = \frac{\mu_0 I}{2\pi \cdot s} \cdot (\sin^2 \alpha_1 + p \cdot \sin^2 \alpha_2)$$

$$B_y = |\vec{B}_1| \cos \alpha_1 + p \cdot |\vec{B}_2| \cos \alpha_2 = \frac{\mu_0 I}{2\pi} \cdot \left(\frac{1}{b+d} \cos^2 \alpha_1 + p \cdot \frac{1}{b} \cdot \cos^2 \alpha_2 \right)$$

For Q on the x-axis, one obtains ($\alpha_1 = \alpha_2 = 0$)

The peak at the minimum of the magnetic field originates from the reflection of the negative magnetic field as positive values, since the measuring instrument only indicates the absolute value of the magnetic field. The different values of the magnetic field at $r = -5 \text{ mm}$ and $r = +5 \text{ mm}$ occur because of the additive or subtractive superimposition of the magnetic fields of conductors 1 and 2.

Results (4/5)




Fig. 6: Magnetic field component B_x of two parallel conductors on the x-axis as a function of the distance from one conductor, if the current in both conductors is in the same direction.




Fig. 7: Magnetic field component B_y of two parallel conductors on the x-axis as a function of the distance from one conductor, if the current in the two conductors is in opposite directions ($I = 107 \text{ A}$).

Results (5/5)

The increase in the field at conductor 2 in comparison with conductor 1 at $r = 65 \text{ mm}$ as compared with $r = 5 \text{ mm}$ occurs because of the higher current density in conductor 2, which results from the resistance of the connecting piece between conductors 1 and 2. Finally, beyond conductor 2 ($r = 75 \text{ mm}$), the effect of conductor 3 becomes noticeable.

This is parallel to conductors 1 and 2, but the current in it flows in the opposite direction to that in conductors 1 and 2 and thus reinforces the magnetic field of 1 and 2 in this area.

The strengthening of the fields can be clearly seen in the space between the two conductors, compared with the reduction in the area beyond the two conductors.