


4 Point Method / Measurement of low resistances / Ohm's Law



Physics

Electricity & Magnetism

Electric current & its effects

Physics

Electricity & Magnetism

Simple circuits, resistors & capacitors

Applied Science

Engineering

Electrical Engineering

Properties of Electrical Circuits



Difficulty level

hard



Group size

2



Preparation time

45+ minutes



Execution time


45+ minutes

This content can also be found online at:



<http://localhost:1337/c/6005692daf145200034371e8>

PHYWE



General information

Application

PHYWE




Fig 1: Experimental set-up

Resistances are very important in electronics and have wide applications in computing and other kind of circuits.

Other information (1/2)

PHYWE

**Prior****knowledge****Main****principle**

The prior knowledge for this experiment is found in the Theory section.

The resistances of various DC conductors are determined by recording the current/voltage characteristic. The resistivity of metal rods and the contact resistance of connecting cords are calculated.

Other information (2/2)

PHYWE

**Learning
objective****Tasks**

The goal of this experiment is to measure low resistances via the 4 point method.

1. To plot the current/voltage characteristics of metal rods (copper and aluminium) and to calculate their resistivity.
2. To determine the resistance of various connecting cords by plotting their current/voltage characteristics and calculating the contact resistances.

Theory (1/2)

PHYWE

The resistivity ρ of the metal is determined from the resistance R of the rod and its dimensions. The rod has a diameter of 2.5 cm (cross section $A = 4.91 \times 10^{-4}$) and is 31.5 cm long (length l) between the two voltmeter connections.

$$\rho = \frac{A \cdot R}{l} \quad (1)$$

Ohm's law

$$U = R \cdot I \quad (2)$$




Fig. 2: Current/voltage characteristics of a copper rod and an aluminium rod.

Theory (2/2)

PHYWE

The regression lines of the measured values in Fig. 2 give

$R_{Cu} = 11.5 \pm 0.3 \mu\Omega$ for the copper rod, and

$R_{Al} = 19.1 \pm 0.2 \mu\Omega$ for the aluminium rod.

The values of resistivity obtained using equation (1) are:

Element measured ρ [$10^{-8} \Omega$] Bibliographic data at 20° ρ [$10^{-8} \Omega$]

Cu	1.79	1.68
Al	2.98	2.72

The aluminium rod is not pure, it contains other additions.

The copper wire in the cords has a cross section A of 2 mm^2

Equipment

Position	Material	Item No.	Quantity
1	Heat conductivity rod, Cu	04518-11	1
2	Heat conductivity rod, Al	04518-12	1
3	PHYWE Universal measuring amplifier	13626-93	1
4	Digital multimeter, 600V AC/DC, 10A AC/DC, 20 M Ω , 200 μ F, 20 kHz, -20°C... 760°C	07122-00	2
5	PHYWE Power supply, 230 V, DC: 0...12 V, 2 A / AC: 6 V, 12 V, 5 A	13506-93	1
6	Connection box	06000-00	1
7	Connecting cord, 100 mm, yellow	07359-02	2
8	Connecting cord, 32 A, 250 mm, red	07360-01	1
9	Connecting cord, 32 A, 250 mm, blue	07360-04	1
10	Connecting cord, 32 A, 500 mm, red	07361-01	2
11	Connecting cord, 32 A, 500 mm, blue	07361-04	1
12	Connecting cord, 32 A, 750 mm, yellow	07362-02	2
13	Connecting cord, 32 A, 750 mm, blue	07362-04	1
14	Connecting cord, 32 A, 2000 mm, yellow	07365-02	2

PHYWE



Setup and Procedure

Setup

PHYWE

1. Connect the metal rod to the mains with an ammeter. Measure the voltage drop across the rod at two sockets on the side, using the amplifier (four-wire method of measurement, see Fig. 1).
2. Connect a connecting cord into the circuit in place of the metal rod, using two double sockets with cross hole (Fig. 3a). Connect the voltmeter to the sockets of the connecting cord connector (similar to the four-wire method; measuring U_1 as shown in Fig. 3). The voltage drops not only across the pure line resistor R_1 but also across the two line/plug contact resistors R_{1p} as well.




Fig. 3a: Measuring the contact resistance and resistivity of connecting cords: sketch of the set-up

Procedure

PHYWE

Determine the total resistance of the connecting cord with connectors by connecting the Voltmeter to the holes in the double sockets (measuring U_2 in Fig. 3). The plug/double socket contact resistances R_{pd} are obtained by comparing U_1 and U_2 .




Fig. 3b: Equivalent circuit diagram: R_t , R_{pd} and R_{1p} are contact resistors, R_1 a line resistor.

PHYWE



Evaluation

Results (1/3)

PHYWE

The line resistance R_1 of the connecting cords can be calculated using (1):

$$R_1 = \rho \cdot \left(\frac{l}{A}\right)$$

The line/plug contact resistance can be established from the difference between the line resistance R_1 calculated and the resistance R_1 measured. R_1 is determined from the slope of the straight lines in Fig. 4.

l [mm]	R_1 [mΩ]	R_1 [mΩ]	$(R_1 - R_1)$ [mΩ]
100	1	0.67	5.6
	2		1.6
750	1	5.0	10.7
	2		9.1
2000	1	13.4	18.6
	2		18.2

Results (2/3)

PHYWE




Fig. 4: Current/voltage characteristics of some connecting cords of different lengths.

The average of the line/plug contact resistance values is:

$$R_{1p} = \frac{R_1 - R_2}{2} = 2.1\text{m}\Omega$$

The plug/double socket contact resistance can be determined by comparing the voltages U_1 and U_2 (see Figs. 3):

$$R_{pd} = \frac{U_1 - U_2}{I}$$

In accordance with Figs. 3b, $U_1 = R_1 \cdot I$ with $R_1 = R_1 + R_{1p_1} + R_{1p_2}$

Results (3/3)

PHYWE

and $U_2 = R_2 \cdot I$ with $R_2 = R_1 + R_{pd_1} + R_{pd_2}$

For a connecting cord 100 mm long the measured values give:

$$R_1 = 5.6m\Omega$$

The plug/double socket contact resistance is therefore of the order of

$$R_{pd} = 30m\Omega$$